

SPRAYING DISINFECTION



ANIOSPRAY QUICK

- » Ethanol-based composition.
- » Aldehyde free.
- » Broad activity spectrum in 30 seconds.
- » Does not leave any trace.
- » Pleasant fragrance.



INSTRUCTIONS FOR USE



READY-TO-USE
Active from 30 sec.

INDICATIONS

Quick disinfection of surfaces, medical equipment and non-immersible and non-invasive medical devices, previously cleaned and alcohol resistant (stethoscopes, cables and connectors, pressure sensors, blood sugar testing material...).

COMPOSITION

Hydroalcoholic solution (ethanol 55%), quaternary ammonium propionate, fragrance.

PACKAGING - REF.

12 bottles, 1L each with sprayers.....2084.073
4 cans, 5L each.....2084.034

ACTIVE AGAINST	ACCORDING TO STANDARDS
BACTERIA	EN 1040, EN 13727, EN 14561, EN 13697 Active against MDRB according to EN 14561, EN 13697
MYCOBACTERIA	EN 14348, EN 14563
YEASTS/MOULDS	EN 1275, EN 13624, EN 14562, EN 13697
VIRUSES	Active against Adenovirus, HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus, Herpesvirus, Coronavirus, Rotavirus, Influenza virus A (H ₁ N ₁), Feline Calicivirus, Norovirus MNV

